

Portable Eyewash Station



EBS 1000

Instructions

- To open pull lightly the transparent lid and it fall downward.
- Take out one of the bottles, twist The Eye Cup cap, the security seal will break enabling you to use the solution.
- The anatomic Eye Cup allows to direct the solution jet correctly, obtaining an efficient eye wash without any waste.
- After washing your eye, the mirror can be useful to find out the eventual damages occurred.



**Spare Bottle 500ml
EEB 500**

- Contains the necessary items for an efficient First Emergency eye wash
- Thermo-formed, has a transparent lid
- Can be easily wall-mounted with screw and plug, supplied in the kit

Contents:

- 2 Sterile Saline Solution of 500 ml
- 2 Integrated Eye cups
- 1 Mirror

Dimensions: 410mm x 320mm x 75mm

Solution Made in Germany

Packed in Italy

EBS 1000
Packing: 1x1set
Dimension: L32 X W9 X H43 CM
Net Weight / KG 1.25 | Gross Weight / KG 1.50

EEB 500
Packing: 1x10nos
Dimension: L30 X W15 X H35 CM
Net Weight / KG 0.55 | Gross Weight / KG 5.50